

> SURGICAL AND EQUIPMENT DRAPES



Code Description Box Qty.

EQUIPMENT DRAPES

SUPREME BACK TABLE COVER

DEF5015	112cm x 228cm (44" x 90")	40
DEF5016	152cm x 228cm (60" x 90")	40
DEF5017	HEAVY DUTY 152cm x 228cm (60" x 90")	40
DEF5018	HEAVY DUTY 178cm x 280cm (70" x 110")	40

MAYO STAND COVERS

DEF5019	SUPREME REINFORCED 127cm x 58.5cm, with Adsorbant pad 50cm x 65cm with Cuff and sealed end	50
DEF761	SEALED ENDS 62cm x 91cm Plastic	100
DEF775C	BLUE PLASTIC WITH CUFF 70cm x 137cm	50



CLEAR PLASTIC DRAPE



DEF742	45cm X 61cm	200
DEF741	61cm X 91cm	25
DEF079	61cm X 91cm	200
DEF749	70cm X 120cm	100
DEF745	76cm X 106cm	100
DEF746	91cm X 150cm	100
DEF752	120cm X 140cm DOUBLE PACKED	100
DEF748	120cm X 140cm	100

Code Description Box Qty.

PATIENT DRAPES

SOFT TOUCH DRAPES

Soft and comfortable, fluid resistant, ideal for minor procedures.

DEF2404	76cm x 106cm no Fenestration	50
---------	------------------------------	----



DEF2412	76cm x 106cm packed in 2's no Fenestration	50
DEF2521	61cm x 61cm packed in 4's no Fenestration	100
DEF2663	120cm x 140cm with 7cm Fenestration in centre of Drape	25
DEF2661	61cm x 61cm with 7cm Circular Fenestration in centre of Drape	50
DEF2681	76cm x 106cm with 7cm Circular Fenestration in centre of Drape	50

POLYBACKED DRAPES

An absorbant cotton top layer, with a waterproof plastic backing. Soft and low linting
Ideal for catheter & minor procedures

DEF1146	40cm x 42cm no Fenestration	25
DEF833	40cm x 42cm no Fenestration	750
DEF2423	40cm x 42cm with 6cm Circular Fenestration	25
DEF1287	40cm x 60cm no Fenestration	100
DEF1288	40cm x 60cm with 10cm Circular Fenestration	100
DEF1289	40cm x 60cm with 10cm x 15cm Fenestration and Adhesive on each corner of drape 2.5cm x 2.5cm	100




DEF1706	50cm x 60cm with 10cm Circular Fenestration	100
DEF2072	75cm x 105cm with 10cm x 15cm Fenestration and 5cm Adhesive length along 10cm sides of Fenestration	25

Code	Description	Box Qty.
------	-------------	----------

AZURE MINOR PROCEDURE DRAPES

- New viscose fabric combined with a PE Laminate to create a drape capable of delivering a high level of performance to most types of surgery.



DEF5090	50cm x 60cm with 6.5 Circular Fenestration with Full Adhesive	25
DEF5091	90cm x 90cm with 10cm x 15cm Fenestration with Full Adhesive	40
DEF5092	200cm x 300cm with 15cm x 25cm Fenestration with Full Adhesive	12

CLEAR PLASTIC EPIDURAL DRAPES WITH 10CM X 15CM FENESTRATION AND 2.5CM ADHESIVE AROUND FENESTRATION

DEF1481	76cm x 106cm	100
---------	--------------	-----



DEF1774	91cm x 115cm	25
DEF748FA-S	120cm x 140cm	10

CLEAR PLASTIC PATIENT DRAPE

DEF2430	91cm x 115cm with 2.5cm x 115cm Adhesive down length of drape	50
---------	---	----

PLASTIC DRAPE YELLOW

DEF381	YELLOW 45mm X 45mm	100
--------	--------------------	-----



Code	Description	Box Qty.
------	-------------	----------

PATIENT SURGICAL DRAPES CLASSIC STOCKINETTE

DEF5060	MEDIUM 20cm x 122cm	25
---------	---------------------	----



DEF5061	LARGE 23cm x 122cm	25
DEF5062	XLARGE 30cm x 122cm	25

SUPREME DRAPE SHEETS

- A Hydrophilic Trilaminar Fabric specifically designed for top level performance coupled with great absorbency, low lint and exceptional patient comfort.
- This fabric is suitable for all levels of surgery.

DEF5024	HALF DRAPE 112cm x 145cm	20
DEF5025	MEDIUM DRAPE 112cm x 180cm	20
DEF5026	3/4 DRAPE 142cm x 193cm	20
DEF5027	3/4 DRAPE 142cm x 193cm with Reinforcement	20
DEF5028	LARGE DRAPE 183cm x 215cm	20
DEF5029	X LARGE DRAPE 178cm x 254cm	12
DEF5030	TOP/HEAD DRAPE 250cm x 135cm with Reinforcement and Adhesive	15
DEF5031	SIDE DRAPE 110cm x 100cm with Reinforcement and Adhesive	20
DEF5032	BOTTOM/FOOT DRAPE 190cm x 185cm with Reinforcement and Adhesive	20
DEF5033	BAR DRAPE 287cm x 183cm with Reinforcement and Adhesive	15

SUPREME IMPERVIOUS SPLITSHEET DRAPE

DEF5049	152cm x 213cm, Split 76cm x 20cm.	25
DEF5050	152cm x 183cm, Split 52cm x 15cm	25

DEF5034	UNIVERSAL PACK	6
---------	-----------------------	---

- 1 x Top Drape 250cm x 135cm. Reinforcement and Adhesive.
- 2 x Side Drapes 110cm x 100cm. Reinforcement and Adhesive.
- 1 x Bottom Drape 190cm x 185cm. Reinforcement and Adhesive
- 1 x Mayo Stand Cover 58cm x 141cm
- 1 x Backtable Cover 152cm x 225cm
- 4 x Utility Drapes

cont...

Code	Description	Box Qty.
DEF5111	UNIVERSAL PACK - SHORT SIDES 1 x Top Drape 250cm x 135cm. Reinforcement and Adhesive. 2 x Side Drape 110cm x 100cm Reinforcement and Adhesive. 1 x Bottom Drape 190cm x 185cm Reinforcement and Adhesive. 1 x Back Table Cover 152cm x 228cm	6
DEF5110	UNIVERSAL PACK - LONG SIDES 1 x Top Drape 250cm x 135cm. Reinforcement and Adhesive. 2 x Side Drapes 110cm x 183cm Reinforcement and Adhesive. 1 x Bottom Drape 190cm x 185cm Reinforcement and Adhesive. 1 x Back Table Cover 152cm x 228cm	6

SPECIALTY PROCEDURE DRAPES**ORTHOPAEDIC DRAPES**

DEF5035	ARTHROSCOPY DRAPE with Armboard Covers and Pouch with Port 447cm x 295cm	6
DEF5036	HIP DRAPE 295cm x 351cm. Bilateral Clear Plastic Leg Pouches Pear Shaped Conformable Fenestration. Armboard Covers. Reinforcement	6
DEF5037	SPLIT U DRAPE 200cm x 300cm. Split 117cm x 15cm (46" x 6") Reinforcement and Adhesive	12
DEF5038	EXTREMITY DRAPE 305cm x 325cm x 227cm. 6cm Conformable Fenestration. Armboard Covers. Reinforcement and Adhesive	6
DEF5039	HAND DRAPE. T Shaped Drape 290cm x 370cm x 190cm. 5cm Circular Conformable Fenestration. Reinforcement and Adhesive.	8
DEF5040	SHOULDER ARTHROSCOPY DRAPE. 255cm x 160cm. Pear Shaped Split Fenestration with Built-In Adhesive Incise. Pouch with Dual Ports.	6
DEF5041	SHOULDER SPLIT DRAPE. 255cm x 160cm. Pear Shaped Split with Adhesive Incise. Reinforcement and Adhesive.	6

GYNAECOLOGY AND OBSTETRIC DRAPES

DEF5042	LSCS DRAPE with built in Incise, 312cm x 249cm. Pear Shaped Fenestration. 360° Pouch: 2 x Ports. Armboard Covers.	8
DEF5043	T SHEET GYNAE DRAPE 140cm x 192cm x 76cm. Perineal Fenestration 8cm x 15cm.	16
DEF5044	LAPAROSCOPY DRAPE 259cm x 246cm x 98cm. Abdominal Triangular Fenestration 19cm x 22cm x 22cm. Perineal Fenestration 8cm x 15cm. Reinforcement.	12
DEF5045	LAPAROSCOPY /LAVH DRAPE 279cm x 276cm x 254cm Abdominal Trapezoidal Fenestration 18cm x 23cm x 23cm. Perineal Fenestration 8cm x 15cm. Attached Leggings, Pouch with Port. Reinforcement. 8/ctn	8
DEF5046	LEGGINGS 73cm x 122cm with 15cm Cuff. Pkt/2	20

Code	Description	Box Qty.
DEF5047	UNDERBUTTOCKS DRAPE 90cm x 120cm with Cuff.	20
DEF5048	UNDERBUTTOCKS DRAPE with Pouch 90cm x 120cm.	20

GENERAL SURGICAL DRAPES

DEF5051	LAPAROTOMY DRAPE with built in Adhesive Incise. 254cm x 314cm x 196cm. Fenestration with Adhesive Incise 23cm x 36cm. Armboard Covers: Reinforcement and Adhesive.	12
DEF5052	ENDOSCOPY DRAPE 254cm x 314cm x 196cm. Fenestration 30cm x 30cm. Two Clear Control Pockets with Velcro Closures. Armboard Covers.	12

CARDIAC AND ANGIOGRAPHY DRAPES

DEF5053	CARDIAC SPLIT 280cm x 350cm x 221cm. Split 17cm x 229cm. Dual Reinforced Troughs with Velcro Fasteners. Armboard Covers. Reinforcement and Adhesive.	6
DEF5054	UNDERLEG PERINEAL DRAPE 213cm x 215cm. Split 30cm x 45cm Reinforcement and Adhesive.	12
DEF5055	ANGIOGRAPHY DRAPE 218cm x 343cm. Dual 9cm Circular Fenestration with Adhesive. Dual Clear Windows. Reinforcement.	10
DEF5056	PACEMAKER DRAPE 222cm x 343cm. Dual 10cm x 15cm. Fenestration with Adhesive Incise. Bilateral Pockets. Clear Breathe Bar Side Windows. Reinforcement.	6

EENT DRAPES

DEF5057	EENT SPLIT DRAPE 196cm x 280cm. Split 6.5cm x 63.5cm. Reinforcement and Adhesive.	12
DEF5058	HEAD TURBAN DRAPE 102cm x 69cm with Adhesive and Turban.	12

OPHTHALMIC DRAPES

DEF5059	140cm x 140cm. 18cm x 15cm Fenestration and Adhesive Incise	8
DEF5063	150cm x 180cm. 6.5cm Circular Absorbant Fenestration	6
DEF5064	EXTRA OCCULAR 160cm x 176cm. 7cm Circular Fenestration	6
DEF5065	166cm x 254cm. Fenestration 10cm x 15.5cm with Adhesive Incise and Pouch	6
DEF5066	186cm x 270cm. Fenestration 9cm x 7cm	6

UROLOGY DRAPES

DEF5067	TUR DRAPE 152cm X 166cm. 6cm Circular Abdominal Aperture. Elastic Perineal Aperture. Non Latex Finger Cot Pouch with Filter and Port. Attached Leggings.	6
---------	--	---