


> SILICONE LOOPS & BOOTS




Code	Description	Box Qty.
------	-------------	----------

IDENTILOOPS

- Silicone radio-opaque Identi Loops are used in surgical situations requiring retraction, identification or occlusion of arteries, veins, nerves, ureters or tendons.
- They are oval in cross-section, allowing for even pressure distribution, minimizing damage to tissues.

DEF1318	2 BLUE, 1 RED, 1 YELLOW, 1 WMAXI-WHITE 40cm 5's	10
---------	--	----



2200620-2	MINI-BLUE 40cm 2's	10
2202620-2	MAXI-BLUE 40cm 2's	10
DEF1282	MAXI-BLUE 40cm 2's D/P	50
DEF1193	SUPER MAXI BLUE 40cm 2's	10



DEF2182	SUPER MAXI BLUE 60cm 2's	10
2200220-2	MINI-RED 40cm 2's	10
2200220-4	MINI-RED 40cm 4's	10
2202220-2	MAXI-RED 40cm 2's	10
2202220-4	MAXI-RED 40cm 4's	10
2200420-2	MINI-YELLOW 40cm 2's	10
2202420-2	MAXI-YELLOW 40cm 2's	10
2200920-2	MINI-WHITE 40cm 2's	10
2202920-2	MAXI-WHITE 40cm 2's	10

Code	Description	Box Qty.
------	-------------	----------

SILCLAMPS

- The Silclamps are covers for surgical clamps, hemostats and forceps.
- They are soft and pliable, providing even pressure distribution, minimizing tissue trauma.

DEF1628	BLUE 3/10" x 10cm 2's	50
DEF1189	RED 1/10" x 10cm SMALL	100

**SUTURE BOOTS SILICONE WITH ADHESIVE TRAY DP**

- Suture Boot Covers offer a non-damaging way to firmly hold delicate sutures without slippage.
- Radio-opaque.

DEF1093	2's	25
---------	-----	----



DEF1094	4's	25
DEF1095	6's	25
DEF1096	8's	25
DEF1097	10's	25
DEF2284	12's	25